

MiniDock MDD

Moleküler Teşhis Cihazı


Teşhiste hız ve güvenilirlik!




Test Prosedürü



Step 1 Numune toplama



Step 2 Karıştırma



Step 3 Kite Ekleme



Step 4 Çalışma

Sonuç Görüntüleme



Mobil uygulama



Işık göstergesi



Bilgisayar yazılımı

Daha fazlası için;



MiniDock MDD

Moleküler Teşhis Cihazı



Pluslife Minidock MDD, küçük ve hafif yapısıyla tanısal testlerde benzersiz bir kullanım kolaylığı ve doğruluk sunan yenilikçi bir moleküler teşhis cihazıdır.

Daha hızlı ve güvenilir sonuçlar için RHAM (Rapid Hybridization Amplification Method) Teknolojisini kullanır.

Özellikleri

- Hızlı ve yüksek doğruluklu sonuçlar
- Sonuç raporu sunabilme
- Telefon ve bilgisayar üzerinden veri aktarımı
- Kullanıcı dostu kullanım
- Geniş test menüsü
- Uygun maliyetli tanı çözümü
- Soğuk zincir gerektirmeyen kitler

Teknik Özellikler

- Boyut : 101x91x65mm
- Ağırlık : 210g
- Prensipl : RHAM Moleküler Teşhis

Parametreler:

FPV/CPV Nucleic Acid Test

CDV/CPIV/CAV-2/Bb Nucleic Acid Test

FeLV/FIV Nucleic Acid Test

Babesia gibsoni/canis Nucleic Acid Test

Anaplasma/Ehrlichia Nucleic Acid Test

Leishmania spp. Nucleic Acid Test

Dirofilaria immitis/Dirofilaria repens Nucleic Acid Test

Giardia lamblia/Tritrichomonas foetus Nucleic Acid Test

FCV/FHV-1/M. felis/C. felis/Bb Nucleic Acid Test

Hemotropic Mycoplasma Mhf/CMhm/CMt Nucleic Acid Test

Giardia/T. foetus/Parvovirus/Coronavirus Nucleic Acid Test